

Guided Acquisition

Advanced Automated Microscopy



Carl Zeiss Microscopy GmbH
Carl-Zeiss-Promenade 10
07745 Jena
Germany

microscopy@zeiss.com
www.zeiss.com/microscopy

Carl Zeiss Microscopy GmbH
ZEISS Group
Kistlerhofstr. 75 81379 München
Germany

Effective from: 01 / 2018

©2018 This document or any part of it must not be translated, reproduced, or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording, or by any information or retrieval system. Violations will be prosecuted. The use of general descriptive names, registered names, trademarks, etc. in this document does not imply, even in the absence of a specific statement, that such names are exempt from the relevant protective laws and regulations and therefore free for general use. Software programs willfully remain the property of ZEISS. No program, documentation, or subsequent upgrade thereof may be disclosed to any third party, unless prior written consent of ZEISS has been procured to do so, nor may be copied or otherwise duplicated, even for the customer's internal needs, apart from a single back-up copy for safety purposes.

ZEISS reserves the right to make modifications to this document without notice.

Contents

1	Introduction	4
1.1	Basic Workflow Definition	5
1.2	General Workflow Definition	6
1.3	General Workflow - Complete Overview	8
2	Workflow - Single Steps	9
2.1	Acquire Sample Data with Overview Scan	9
2.2	Create Image Analysis Pipeline	10
2.3	Define Overview Experiment	13
2.4	Define Detailed Scan Experiment	14
3	Automation via OAD	15
3.1	Prerequisites	15
3.2	Pure Scripting or User Dialog	15
3.3	User Dialog Requirements	15
3.4	Prerequisites	16
3.5	Create the User Dialog	19
3.6	Overview Scan Experiment	22
3.7	Find the interesting Objects	23
3.8	Modify the Tile Dimensions	25
3.9	Run the Detailed Scan	26
4	Testrun with FluoCells Slide	29
5	Test Run with Brain Slide	31
6	Mitosis Detection with Camera and LSM	33
7	Appendix	35
7.1	Complete Workflow Diagram	35
7.2	Python Script: Guided_Acquisition_shortUI.py	36
8	ToDo's and Limitations	43
9	Disclaimer	43

1 Introduction

For a growing number of applications, it will be crucial to acquire data in a smart way. One way to achieve this goal is to build a smart microscope, which essentially means creating smart software workflows to control the hardware based on image analysis results.

Idea or Task:

- Scan or inspect a large area.
- Detect an "interesting" object.
- Acquire detailed data for every event.
- Automate the workflow to minimize user interference.

First of all it is important to define what a **Object of Interest** can actually be. Example would be an object that meets specific criteria, for example:

- size
- brightness
- shape
- intensity
- combinations of the above

It could be something quite simple. For instance one can have lots of cells, that are stained with blue dye, and only a few of them (maybe where the transfection worked ...) are also expressing GFP. The idea here would be to detect all cells that are positive for both colors and acquire an z-Stack for every cell (blue & green) that meets those criteria. Therefore this kind of application requires three major tasks:

1. **Define the Overview Scan Experiment.**
2. **Define the object detection rules, e.g. setup image analysis.**
3. **Define the Detailed Scan(s) to be carried out in case of a "positive" object.**

1.1 Basic Workflow Definition

The main workflow can be summarized as follows:

1. Acquire some sample data showing a object of interest.
2. Setup an Image Analysis Pipeline to detect those events.
3. Define an experiment which does the Overview Scan.
4. Define an experiment which does the Detailed Scan.
5. Automate the entire workflow.

The goal of this tutorial is to create an automated workflow that can be used to easily setup a **Guided Acquisition**. This requires some knowledge about the OAD macro environment and its scripting language Python.

1.2 General Workflow Definition

This is the outline of the general workflow, which would be automated.




Figure 1: General workflow for Guided Acquisition.

A different way to visualize this is shown in the figure below.




Figure 2: Larger area with some candidates for a "interesting" objects.

1.3 General Workflow - Complete Overview




Figure 3: Workflow - Actionable Information 1.




Figure 4: Workflow - Actionable Information 2

2 Workflow - Single Steps

This part of the tutorial will explain all required steps to set up the complete workflow in more detail. Some knowledge about image analysis is required here.

2.1 Acquire Sample Data with Overview Scan

The first crucial step is to have a basic idea of what the actual rare event will look like, inside an real image acquired, with the appropriate parameters (light intensity, detector settings, filters, objective, ...). An ideal sample data set should be acquired using the "real" acquisition parameters will be used to define the overview scan experiment later on.

Typically such an overview scan is acquired using a lower magnification in combination with a tile experiment. The exact parameters will be specific for the application, but the main idea is always the same:

The overview image must contain the information to locate the objects of interest based on "some" features that can be retrieved via an appropriate image analysis.

Once a representative sample data set is available, one can start setting up an image analysis pipeline. ZEN offers currently two ways do to so:

1. The easy way – Use the Analysis Wizard.
2. Program your own image analysis pipeline using an OAD macro (this is not covered by this tutorial).
3. Program your own image analysis using an external software (e.g. Fiji, Python etc.) and use it from within ZEN (this is not covered by this tutorial).

The most comfortable way is to use the wizard; this offers enough flexibility to cover most of the applications.

2.2 Create Image Analysis Pipeline

Here one already sees a "half-ready" image analysis pipeline. In order to automate, the wizard step containing the feature selection is crucial.

In order to relocate all detected objects, it is of course required to retrieve the current XY(Z) position of every detected event.

The main idea is to determine all parameters required to get access to the stage coordinates of a detected object.




Figure 5: Select the required positional features inside the wizard.

Currently one needs to select the following parameters:

- **Object ID**
- **Bound Center Stage X**
- **Bound Center Stage Y**
- **Bound Width**
- **Bound Height**
- **ImageSceneContainerName**

The object ID is needed to keep track of the executed detail scans and to save the resulting image files with the correct ID as part of the used filename.

The parameters **Bound Center Stage X** and **Bound Stage Center Y** yield absolute stage coordinates of the detected objects. These will be used to relocate the objects for the detailed scan later on.

The parameters **Bound Width** and **Bound Height** yield in the absolute width and height of the bounding rectangle. This is required if the detected object is bigger than a single frame. In such a case the tile region of the detailed experiment will be adapted to the size of the bounding rectangle automatically.

The parameter **ImageSceneContainerName** is needed to keep track of the respective wells, because the wellID will be used as part of the filename for the detail scan. Depending on your application, it might be useful to measure additional parameters. Feel free to add whatever might be useful.

INFO:

The parameter **ImageSceneContainerName** is not required under all circumstances and only makes sense if there really is a physical well or container on your sample.

2.3 Define Overview Experiment

Usually this is done using the Tiles and Positions module to scan a large area. This tutorial assumes that one is already familiar with this ZEN module.

Important Remark: Since one wants to relocate the detected objects, it is required to calibrate the XY stage before the overview scan, and a rigid and stiff sample holder should be used.

For this tutorial one can use an already acquired image representing the real overview scan (instead of a real experiment). The shown image is a 16x13 tile image acquired with a 10X magnification. The result of such a scan followed by the image analysis is shown here:




Figure 6: Result of the image analysis on the overview scan image.

2.4 Define Detailed Scan Experiment

First of all it is important to understand, that there is no such thing as a typical "Detailed Scan". What this experiment (or even a workflow) might be, depends on the application. In a general sense such a detailed scan is "some kind" of experiment carried at a specific position based upon the results of the image analysis. Examples could be:

- Simple Z-Stack with a high-NA objective lens.
- Multi-Channel Z-Stack using an optical sectioning method like SD, Apotome or AiryScan
- One of the above combined with a tile experiment.
- A series of subsequent experiments with different image modalities.

Since this concept is based on OAD, it is possible to automate almost any kind of workflow at a given position.




Figure 7: Setup of a simple detailed scan experiment.

Now the preparation is complete. So far we have available:

- The image analysis pipeline based on some sample data
- The Overview Scan Experiment = tile experiment created by the Tiles & Positions Module
- The Detailed Scan Experiment, which is also a tile experiment.

And now the real interesting part begins – how to setup a complete workflow out of those building blocks using an OAD Python script inside ZEN Blue.

3 Automation via OAD

3.1 Prerequisites

In order to be able to run the script one needs to copy it into the correct folder, which is located at:

"c:/Users/XXXX/Documents/Carl Zeiss/ZEN/Documents/Macros"

where XXXX stand for the actual user name.

3.2 Pure Scripting or User Dialog

First of all one has to decide if a pure script is sufficient or if there should be some kind of simple user interface (very basic wizard).

- A pure script is of course the faster solution, but is may be difficult to use for less advanced users.
- A basic wizard offers less flexibility, but may be easier to use for most users.

The choice clearly depends on the user, since there is no right or wrong approach here. For this tutorial we will focus on the 2nd approach. Once one has figured out how to realize this, the 1st approach will be straight forward since the basic workflow is identical.

3.3 User Dialog Requirements

The basic user interface of the dialog should have the following functionality:

- **Select the overview scan experiment.**
- **Select the image analysis pipeline** (to analyze the overview image)
- **Select the detailed scan experiment** (to acquire at "found" positions)
- **Define options for focussing** (depends on available hardware)
- **Select a folder to store the image data and analysis results.**

3.4 Prerequisites

Now we are ready to create the OAD macro required to control the complete workflow. The first step is (as always) to import the required modules and define the folder locations, etc. Additionally there are several python functions defined here for later usage.

```
1 ######
2 # File      : Guided_Acquisition_shortUI.py
3 # Version   : 7.3
4 # Author    : czsrh, czmla, czkel
5 # Date      : 19.08.2019
6 # Institution : Carl Zeiss Microscopy GmbH
7 #
8 # !!! Requires with ZEN >=2.6 HF3 - Use at your own Risk !!!
9 #
10 # Optimized for the use with Celldiscoverer 7 and DF2, but
11 # applicable for all motorized stands running in ZEN Blue.
12 # Please adapt focussing commands, especially FindSurface
13 # when using with other stands.
14 #
15 # 1) - Select Overview Scan Experiment
16 # 2) - Select appropriate Image Analysis Pipeline
17 # 3) - Select Detailed Scan Experiment
18 # 4) - Specify the output folder for the image and data tables
19 #
20 #
21 # Copyright(c) 2019 Carl Zeiss AG, Germany. All Rights Reserved.
22 #
23 # Permission is granted to use, modify and distribute this code,
24 # as long as this copyright notice remains part of the code.
25 #####
26
27 import time
28 from datetime import datetime
29 import errno
30 from System import Array
31 from System import ApplicationException
32 from System import TimeoutException
33 from System.IO import File, Directory, Path
34 import sys
35
36
37 # version number for dialog window
38 version = 7.3
39 # file name for overview scan
40 ovscan_name = 'OverviewScan.czii'
41
42 """
43 Additional XY offset for possible 2nd port relative to the 1st port
44 Example: 1st port Camera and 2nd port LSM
45 Can be set to Zero, when system is correctly calibrated
46 !!! Only use when overview and detailed scan are using different detector !!!
47 """
48 dx_detector = 0.0
49 dy_detector = 0.0
50
51 # experiment blockindex
52 blockindex = 0
53 # delay for specific hardware movements in [seconds]
54 hwdelay = 1
55 # posprocessing switch
56 do_postprocess = False
57
58
59 def dircheck(basefolder):
60
61     # check if the destination basefolder exists
62     base_exists = Directory.Exists(basefolder)
63
64     if base_exists:
65         print('Selected Directory Exists: ', base_exists)
66         # specify the desired output format for the folder, e.g. 2017-08-08_17-47-41
67         format = '%Y-%m-%d_%H-%M-%S'
68         # create the new directory
69         newdir = createfolder(basefolder, formatstring=format)
70         print('Created new directory: ', newdir)
71     if not base_exists:
72         Directory.CreateDirectory(basefolder)
73         newdir = basefolder
74
75 return newdir
76
```

```

77 def createfolder(basedir, formatstring='%Y-%m-%d_%H-%M-%S'):
78     # construct new directory name based on date and time
79     newdir = Path.Combine(basedir, datetime.now().strftime(formatstring))
80     # check if the new directory (for whatever reasons) already exists
81     try:
82         newdir_exists = Directory.Exists(newdir)
83         if not newdir_exists:
84             # create new directory if it does not exist
85             Directory.CreateDirectory(newdir)
86         if newdir_exists:
87             # raise error if it really already exists
88             raise SystemExit
89     except OSError as e:
90         if e.errno != errno.EEXIST:
91             newdir = None
92             raise # This was not a "directory exist" error..
93
94     return newdir
95
96
97 def getshortfiles(filelist):
98     files_short = []
99     for short in filelist:
100        files_short.append(Path.GetFileName(short))
101
102     return files_short
103
104
105 def cloneexp(expname, prefix='GA_', save=True, reloadexp=True):
106
107     exp = Zen.Acquisition.Experiments.GetByName(expname)
108     exp_newname = prefix + expname
109
110     # save experiment
111     if save:
112         exp.SaveAs(exp_newname, False)
113         print('Saved Temporary Experiment as : ', exp_newname)
114         # close the original experiment object
115         exp.Close()
116         time.sleep(1)
117
118     # reload experiment
119     if reloadexp:
120         exp_reload = Zen.Acquisition.Experiments.GetByName(exp_newname)
121     elif not reloadexp:
122         exp_reload = None
123
124     return exp_reload
125
126
127 def runSWAF_special(SWAF_exp,
128                      delay=5,
129                      searchStrategy='Full',
130                      sampling=ZenSoftwareAutofocusSampling.Coarse,
131                      relativeRangeIsAutomatic=False,
132                      relativeRangeSize=500,
133                      timeout=0):
134
135     # get current z-Position
136     zSWAF = Zen.Devices.Focus.ActualPosition
137     print('Z-Position before special SWAF : ', zSWAF)
138
139     # set DetailScan active and wait for moving hardware due to settings
140     SWAF_exp.SetActive()
141     time.sleep(delay)
142
143     # set SWAF parameters
144     SWAF_exp.SetAutofocusParameters(searchStrategy=searchStrategy,
145                                     sampling=sampling,
146                                     relativeRangeIsAutomatic=relativeRangeIsAutomatic,
147                                     relativeRangeSize=relativeRangeSize)
148
149     try:
150         print('Running special SWAF ...')
151         zSWAF = Zen.Acquisition.FindAutofocus(SWAF_exp, timeoutSeconds=timeout)
152     except ApplicationException as e:
153         print('Application Exception : ', e.Message)
154     except TimeoutException as e:
155         print(e.Message)
156
157     print('Z-Position after initial SWAF : ', zSWAF)
158
159     return zSWAF
160
161 def checktableentry(datatable, entry2check='ImageSceneContainerName'):
162
163     num_col = datatable.ColumnCount
164     entry_exists = False
165     column = None

```

```
166
167     for c in range(0, num_col):
168         # get the current column name
169         colname = datatable.Columns[c].ColumnName
170         if colname == entry2check:
171             column = c
172             entry_exists = True
173             break
174
175     return entry_exists, column
```

3.5 Create the User Dialog

The next task is to create the user dialog using the requirements referred to earlier on. Once the dialog is closed, the results have to be collected.

Additionally one must check if the Detailed Scan experiment is a tile experiment, where the tile size has to be adapted accordingly. This is done using a little tool function.

```
217 # Initialize Dialog
218 GuidedAcqDialog = ZenWindow()
219 GuidedAcqDialog.Initialize('Guided Acquisition - Version : ' + str(version))
220 # add components to dialog
221 GuidedAcqDialog.AddLabel('1) Select Overview Experiment -----')
222 GuidedAcqDialog.AddDropDown('overview_exp', 'Overview Scan Experiment', expfiles_short, 0)
223 GuidedAcqDialog.AddCheckbox('fs_before_overview', 'OPTION - FindSurface (DF only) before Overview', False)
224 GuidedAcqDialog.AddCheckbox('SWAF_before_overview', 'OPTION - SWAF before Overview', False)
225 GuidedAcqDialog.AddIntegerRange('SWAF_ov_initial_range', 'Initial SWAF Range before Overview [micron]', 200, 50, 3000)
226 GuidedAcqDialog.AddLabel('2) Select Image Analysis to detect objects -----')
227 GuidedAcqDialog.AddDropDown('ip_pipe', 'Image Analysis Pipeline', ipfiles_short, 0)
228 GuidedAcqDialog.AddLabel('3) Select DetailScan Experiment -----')
229 GuidedAcqDialog.AddDropDown('detailed_exp', 'Detailed Scan Experiment', expfiles_short, 1)
230 GuidedAcqDialog.AddCheckbox('fs_before_detail', 'OPTION - FindSurface (DF only) before Detail', False)
231 GuidedAcqDialog.AddCheckbox('SWAF_before_detail', 'OPTION - SWAF before Detail', False)
232 GuidedAcqDialog.AddIntegerRange('SWAF_detail_initial_range', 'Initial SWAF Range before Detail [micron]', 100, 10, 1000)
233 GuidedAcqDialog.AddCheckbox('recallfocus_beforeDT', 'OPTION - Use RecallFocus (DF only) before Detail', False)
234 GuidedAcqDialog.AddLabel('4) Specify basefolder to save the images -----')
235 GuidedAcqDialog.AddFolderBrowser('outfolder', 'Basefolder for Images and Data Tables', imgfolder)
236
237 # show the window
238 result = GuidedAcqDialog.Show()
239 if result.HasCanceled:
240     message = 'Macro was canceled by user.'
241     print(message)
242     raise SystemExit
243
```

The resulting dialog window is shown below.




Figure 8: Simple User Dialog to run the Guided Acquisition WorkFlow (Remark: Version in screenshot may differ from the actual latest version)

The dialog has four main steps where the user must define the desired options:

1. Select **overview experiment** and additional focus options.
2. Select **image analysis** to detect the objects.
3. Select **detailed experiment** and additional focus options.

4. Choose an **output folder** to save the data.

3.6 Overview Scan Experiment

At this point we have all the information we need to actually start the workflow. The first step is to execute the actual overview scan experiment, which is most likely a tile experiment. The resulting image will be saved to the specified folder.

```
277 ##### START OVERVIEW SCAN EXPERIMENT #####
278
279 if fs_beforeOV:
280     # initial focussing via FindSurface to assure a good starting position
281     Zen.Acquisition.FindSurface()
282     print('Z-Position after FindSurface: ', Zen.Devices.Focus.ActualPosition)
283
284 if SWAF_beforeOV:
285     zSWAF = runSWAF_special(OVScan_reloaded,
286                             delay=hwdelay,
287                             searchStrategy='Full',
288                             sampling=ZenSoftwareAutofocusSampling.Coarse,
289                             relativeRangeIsAutomatic=False,
290                             relativeRangeSize=SWAF_beforeOV_range,
291                             timeout=1)
292
293 # get the resulting z-position
294 znew = Zen.Devices.Focus.ActualPosition
295
296 # adapt the Overview Scan Tile Experiment with new Z-Position
297 if OVScanIsTileExp:
298     OVScan_reloaded.ModifyTileRegionsZ(blockindex, znew)
299     print('Adapted Z-Position of Tile OverView. New Z = ', '%.2f' % znew)
300
301 # execute the experiment
302 print('\nRunning OverviewScan Experiment.\n')
303 output_OVScan = Zen.Acquisition.Execute(OVScan_reloaded)
304 # For testing purposes - Load overview scan image automatically instead of executing the "real" experiment
305 # output_OVScan = Zen.Application.LoadImage(r'c:\Temp\input\OverViewScan_8Brains.czi', False)
306
307 # show the overview scan inside the document area
308 Zen.Application.Documents.Add(output_OVScan)
309 ovdoc = Zen.Application.Documents.GetByName(output_OVScan.Name)
310
311 # save the overview scan image inside the select folder
312 output_OVScan.Save(Path.Combine(OutputFolder, ovscan_name))
313
314 ##### END OVERVIEW SCAN EXPERIMENT #####
315
```

A very important part of such an automated workflow is to find the focus. To ensure this, the desired objective and optovar are changed before the overview starts, followed by a series of **optional** focus actions:

1. Find the surface of the sample carrier.
2. Focus onto the specimen sample via SWAF.
3. Store the resulting new Z-position.

When testing out applications like this it is very useful to use test data sets to save time instead of re-running an acquisition many times. Therefore it can be helpful to **comment** the lines where the real experiment is executed and **uncomment** the line inside the script, where one can load an already acquired overview image.

3.7 Find the interesting Objects

Now it is time to run the image analysis on the overview image. As a result, two tables will be created inside ZEN. The first contains information about all objects and the second contains information about the single objects (for instance their stage positions). This second table is what will be used as an input for the detailed scan later on.

```
315 # Load analysis setting created by the wizard or an separate macro
316 ias = ZenImageAnalysisSetting()
317
318 # for simulation use: 000 - RareEventExample.czias
319 ias.Load(ImageAS)
320
321 # Analyse the image
322 Zen.Analyzing.Analyze(output_OVScan, ias)
323
324 # Create Zen table with results for all detected objects (parent class)
325 AllObj = Zen.Analyzing.CreateRegionsTable(output_OVScan)
326
327 # Create Zen table with results for each single object
328 SingleObj = Zen.Analyzing.CreateRegionTable(output_OVScan)
329
330 # check for existence of required column names inside table
331 soi = SingleObj.GetColumnInfoFromImageAnalysis(True)
332
333 # 1st item is a bool indicating if all required columns could be found
334 columnsOK = soi.AreRequiredColumnsAvailable
335
336
337 if not columnsOK:
338     print('Execution stopped. Required Columns are missing.')
339     raise Exception('Execution stopped. Required Columns are missing.')
340
341 # show and save data tables to the specified folder
342 Zen.Application.Documents.Add(AllObj)
343 Zen.Application.Documents.Add(SingleObj)
344 AllObj.Save(Path.Combine(OutputFolder, 'OverviewTable.csv'))
345 SingleObj.Save(Path.Combine(OutputFolder, 'SingleObjectsTable.csv'))
346
347 # check the number of detected objects = rows inside image analysis table
348 num_POI = SingleObj.RowCount
```

As shown in figure 5, it is required to have the correct image analysis features selected. Only the will the respective columns inside the table exist. This has to be checked because otherwise the workflow will be impossible.

The "checking" itself is again done by a tool function. The output will be a boolean and a dictionary in order to look up the column index for the respective image analysis parameters. This has the advantage of being independent from a specific column order. One just has to make sure that all required columns are there.

If the option inside the initial dialog window was checked, the offset between the sample carrier surface and the actual specimen will be measured and stored inside the Definite Focus via **Store Focus**, so that **Recall Focus** can be used later on when needed. The dialog also offers the possibility to enter an arbitrary offset manually, which will only be applied if the offset was not determined automatically.




Figure 9: The results of the overview scan and the image analysis.

Finally the number of detected objects is checked by looking up the number of rows inside the table containing the information about the single objects.

3.8 Modify the Tile Dimensions

Once one has the list of objects with the stage coordinates, the next task is to create a loop to cycle through all detected positions.

- Move to the correct XYZ stage positions.
- Run Focus Actions when required.
- Initialize the DetailScan experiment at the new position.
- Adapt TilePosition using Tiles.TileTools.ModifyTileRegionsXYZ.
- Modify TileSize using Tiles.TileTools.ModifyTileRegionsSize.

```

406 ##### START DETAILED SCAN EXPERIMENT #####
407
408 # get the actual Focus position
409 zpos = Zen.Devices.Focus.ActualPosition
410
411 # check for the column 'ID' which is required
412 ID_exists, column_ID = checktableentry(SingleObj,
413                                         entry2check='ID')
414
415
416 # execute detailed experiment at the position of every detected object
417 for i in range(0, num_POI, 1):
418
419     # get the object information from the position table
420     # get the ID of the object - IDs start with 2 !!!
421     #POI_ID = SingleObj.GetValue(i, 0)
422     POI_ID = SingleObj.GetValue(i, column_ID)
423
424     # get XY-stage position from table
425     xpos = SingleObj.GetValue(i, soi.CenterXColumnIndex)
426     ypos = SingleObj.GetValue(i, soi.CenterYColumnIndex)
427
428     # move to the current position
429     Zen.Devices.Stage.MoveTo(xpos + dx_detector, ypos + dy_detector)
430     print('Moving Stage to Object ID:', POI_ID, ' at : ', '%.2f' % xpos, '%.2f' % ypos)
431
432     # try to apply RecallFocus (DF only) when this option is used
433     if userecallfocus:
434         try:
435             # apply RecallFocus for the current position
436             Zen.Acquisition.RecallFocus()
437             zpos = Zen.Devices.Focus.ActualPosition
438             print('Recall Focus (Definite Focus) applied.')
439             print('Updated Z-Position: ', zpos)
440         except ApplicationException as e:
441             print('Application Exception : ', e.Message)
442             print('RecallFocus (Definite Focus) failed.')
443
444         print('New Z-Position before Detail Experiment will start:', zpos)
445
446     # if DetailScan is a Tile Experiment
447     if DetailIsTileExp:
448
449         print('Detailed Experiment contains TileRegions.')
450         # Modify tile center position - get bounding rectangle width and height in microns
451         bcwidth = SingleObj.GetValue(i, soi.WidthColumnIndex)
452         bcheight = SingleObj.GetValue(i, soi.HeightColumnIndex)
453         print('Width and Height : ', '%.2f' % bcwidth, '%.2f' % bcheight)
454         print('Modifying Tile Properties XYZ Position and width and height.')
455
456         # Modify the XYZ position and size of the TileRegion on-the-fly
457         print('Starting Z-Position for current Object: ', '%.2f' % zpos)
458         print('New Tile Properties: ', '%.2f' % xpos, '%.2f' % ypos, '%.2f' % zpos, '%.2f' % bcwidth, '%.2f' % bcheight)
459         DetailScan_reloaded.ClearTileRegionsAndPositions(blockindex)
460         try:
461             DetailScan_reloaded.AddRectangleTileRegion(blockindex, xpos, ypos, bcwidth, bcheight, zpos)
462         except ApplicationException as e:
463             print('Application Exception : ', e.Message)
464
465     if not DetailIsTileExp:
466         print('Detailed Experiment does not contain TileRegions. Nothing to modify.')
467

```

3.9 Run the Detailed Scan

The last steps, which need to be executed now, are:

- Run the (modified) detailed scan experiment here.
- Save the data to the specified folder.
- Close the image document in ZEN.
- Rename the file using the current object ID.

```
468 # execute the DetailScan experiment
469 print('Running Detail Scan Experiment at new XYZ position.')
470 try:
471     output_detailscan = Zen.Acquisition.Execute(DetailScan_reloaded)
472 except ApplicationException as e:
473     print('Application Exception : ', e.Message)
474
475 # get the image data name
476 dtscan_name = output_detailscan.Name
477
478 """
479 Modification for multiscene images: container name needs to be part
480 of the filename of the detailed Scan.
481 First check if Column "Image Scene Container Name" is available
482 and then add Container Name to filename.
483 """
484
485 wellid_exist, column_wellid = checktableentry(SingleObj,
486                                             entry2check='ImageSceneContainerName')
487
488 # in case Image Scene Container Name is not defined
489 if not wellid_exist:
490     message = 'Missing column ImageSceneContainerName in Image Analysis Results.\nPlease Select Features inside the Image
491     ↪ Analysis when needed.'
492     print(message)
493     container_name = 'empty'
494
495 # save the image data to the selected folder and close the image
496 output_detailscan.Save(Path.Combine(OutputFolder, output_detailscan.Name))
497 output_detailscan.Close()
498
499 # rename the CZI regarding to the object ID - Attention - IDs start with 2 !!!
500 if not wellid_exist:
501     # no wellID found inside table
502     newname_dtscan = 'DTScan_ID_' + str(POI_ID) + '.czi'
503 if wellid_exist:
504     # wellID was found inside table
505     well_id = SingleObj.GetValue(i, column_wellid)
506     newname_dtscan = 'DTScan_Well_' + str(well_id) + '_ID_' + str(POI_ID) + '.czi'
507
508 print('Renaming File: ' + dtscan_name + ' to: ' + newname_dtscan + '\n')
File.Move(Path.Combine(OutputFolder, dtscan_name), Path.Combine(OutputFolder, newname_dtscan))
```




Figure 10: The results of the detailed scan - tile regions adapts to objects.

That is it. All image data sets acquired at the detected positions are stored inside the correct folder using the object IDs as names. The tables are saved in the same location.




Figure 11: The results of the detailed scan stored inside a folder.

4 Testrun with FluоСells Slide

To test the workflow on a real sample one can use the FluоСells slide. The general idea is to setup an image analysis that detects some "interesting" objects. In order to create a somewhat short test run it is recommended to use an image analysis pipeline that only yields in a few positive detection events.




Figure 12: The results of the overview scan and image analysis.




Figure 13: Detailed view on a detected object.




Figure 14: Detailed images from some of the detected objects.

The image analysis yielded **24** positive detections of cells with a large nucleus, preferably cells which were "caught" during cell division. Keep in mind, that the image analysis was not "fine-tuned" to catch all cell divisions. It is all about testing out the general workflow. For this reason it is also good practice to acquire a small overview scan at the beginning, e.g. a 2x2 tile image.

5 Test Run with Brain Slide

The same basic workflow can be also used to detect brain slices on a slide. All one needs to change is the image analysis pipeline in order to detect the slices.




Figure 15: The results of the overview scan and image analysis.




Figure 16: A more detailed image of the first specimen with three zoomed-in regions.

The image analysis yielded **eight** positive regions for brains slices on this slide. The properties of the bounding rectangle were used to modify the required tile experiment. The final detailed scan contained 36 tiles and it was imaged using only a hardware-based focus system.

Keep in mind that the actual focusing during such a experiment must be still part of the focus strategies with the Detail Scan experiment. Which one must be used here greatly depends on the nature of the sample and the actual application. The Guided Acquisition tool only supports finding the initial positions based on the overview scan image.

6 Mitosis Detection with Camera and LSM

An especially interesting option is to combine the power of a camera-based overview scan with the optional sectioning capabilities of an LSM. Such a workflow can be easily configured in ZEN Blue by setting up the respective experiment for the overview scan with camera detection using a low magnification and the LSM-based detail scan, e.g. Z-Stack, using a high NA objective.

Since ZEN Blue 2.5 this software has a module called **ZEN Connect**, which allows combining and correlating images inside one sample-centric workspace. Every acquired image by either the camera or the LSM will be placed here based on the XY stage coordinates. Therefore the **Correlative Workspace (CWS)** is ideally suited to display the results of a Guided Acquisition workflow.




Figure 17: The results of the overview scan and image analysis inside the CWS.

Advanced Automated Microscopy: Guided Acquisition




Figure 18: Overlay between positive detection inside overview scan and detailed scan using the Correlative Workspace (CWS).




Figure 19: Side-by-Side view of the CWS and the normal view inside ZEN Blue.

7 Appendix

7.1 Complete Workflow Diagram

This is the complete workflow diagram for **Guided Acquisition**.




Figure 20: Complete workflow diagram for Guided Acquisition

7.2 Python Script: Guided_Acquisition_shortUI.py

This is the complete ZEN Blue python script used to realize the workflow described above. In principle it can work on any motorized microscope stand running under ZEN blue (with some modifications), but it was tested and optimized mainly for the Celldiscoverer 7 and for the AxioObserver 7. Both of those stands have the Definite Focus 2 as an additional hardware option, which is strongly advised for such highly automated workflows, where fast and efficient focus options, such as Find Surface are really a big plus.

```
1 #####  
2 # File      : Guided_Acquisition_shortUI.py  
3 # Version   : 7.3  
4 # Author    : czsrh, czmla, czkel  
5 # Date      : 19.08.2019  
6 # Institution : Carl Zeiss Microscopy GmbH  
7 #  
8 # !!! Requires with ZEN >=2.6 HF3 - Use at your own Risk !!!  
9 #  
10 # Optimized for the use with Celldiscoverer 7 and DF2, but  
11 # applicable for all motorized stands running in ZEN Blue.  
12 # Please adapt focussing commands, especially FindSurface  
13 # when using with other stands.  
14 #  
15 # 1) - Select Overview Scan Experiment  
16 # 2) - Select appropriate Image Analysis Pipeline  
17 # 3) - Select Detailed Scan Experiment  
18 # 4) - Specify the output folder for the image and data tables  
19 #  
20 #  
21 # Copyright(c) 2019 Carl Zeiss AG, Germany. All Rights Reserved.  
22 #  
23 # Permission is granted to use, modify and distribute this code,  
24 # as long as this copyright notice remains part of the code.  
25 #####  
26  
27 import time  
28 from datetime import datetime  
29 import errno  
30 from System import Array  
31 from System import ApplicationException  
32 from System import TimeoutException  
33 from System.IO import File, Directory, Path  
34 import sys  
35  
36  
37 # version number for dialog window  
38 version = 7.3  
39 # file name for overview scan  
40 ovscan_name = 'OverviewScan.czi'  
41  
42 """  
43 Additional XY offset for possible 2nd port relative to the 1st port  
44 Example: 1st port Camera and 2nd port LSM  
45 Can be set to Zero, when system is correctly calibrated  
46 !!! Only use when overview and detailed scan are using different detector !!!  
47 """  
48 dx_detector = 0.0  
49 dy_detector = 0.0  
50  
51 # experiment blockindex  
52 blockindex = 0  
53 # delay for specific hardware movements in [seconds]  
54 hwdelay = 1  
55 # postprocessing switch  
do_postprocess = False  
56  
57  
58 def dircheck(basefolder):  
59  
60     # check if the destination basefolder exists  
61     base_exists = Directory.Exists(basefolder)  
62  
63     if base_exists:  
64         print('Selected Directory Exists: ', base_exists)  
65         # specify the desired output format for the folder, e.g. 2017-08-08_17-47-41  
66         format = '%Y-%m-%d_%H-%M-%S'  
67         # create the new directory  
68         newdir = createfolder(basefolder, formatstring=format)  
69         print('Created new directory: ', newdir)
```

```

71     if not base_exists:
72         Directory.CreateDirectory(basefolder)
73         newdir = basefolder
74
75     return newdir
76
77
78 def createfolder(basedir, formatstring='%Y-%m-%d %H-%M-%S'):
79     # construct new directory name based on date and time
80     newdir = Path.Combine(basedir, datetime.now().strftime(formatstring))
81     # check if the new directory (for whatever reasons) already exists
82     try:
83         newdir_exists = Directory.Exists(newdir)
84         if not newdir_exists:
85             # create new directory if it does not exist
86             Directory.CreateDirectory(newdir)
87         if newdir_exists:
88             # raise error if it really already exists
89             raise SystemExit
90     except OSError as e:
91         if e.errno != errno.EEXIST:
92             newdir = None
93             raise # This was not a "directory exist" error..
94
95     return newdir
96
97
98 def getshortfiles(filelist):
99     files_short = []
100    for short in filelist:
101        files_short.append(Path.GetFileName(short))
102
103    return files_short
104
105
106 def cloneexp(expname, prefix='GA_', save=True, reloadexp=True):
107
108     exp = Zen.Acquisition.Experiments.GetByName(expname)
109     exp_newname = prefix + expname
110
111     # save experiment
112     if save:
113         exp.SaveAs(exp_newname, False)
114         print('Saved Temporary Experiment as : ', exp_newname)
115         # close the original experiment object
116         exp.Close()
117         time.sleep(1)
118     # reload experiment
119     if reloadexp:
120         exp_reload = Zen.Acquisition.Experiments.GetByName(exp_newname)
121     elif not reloadexp:
122         exp_reload = None
123
124     return exp_reload
125
126
127 def runSWAF_special(SWAF_exp,
128                      delay=5,
129                      searchStrategy='Full',
130                      sampling=ZenSoftwareAutofocusSampling.Coarse,
131                      relativeRangeIsAutomatic=False,
132                      relativeRangeSize=500,
133                      timeout=0):
134
135     # get current z-Position
136     zSWAF = Zen.Devices.Focus.ActualPosition
137     print('Z-Position before special SWAF : ', zSWAF)
138
139     # set DetailScan active and wait for moving hardware due to settings
140     SWAF_exp.SetActive()
141     time.sleep(delay)
142
143     # set SWAF parameters
144     SWAF_exp.SetAutofocusParameters(searchStrategy=searchStrategy,
145                                     sampling=sampling,
146                                     relativeRangeIsAutomatic=relativeRangeIsAutomatic,
147                                     relativeRangeSize=relativeRangeSize)
148
149     try:
150         print('Running special SWAF ...')
151         zSWAF = Zen.Acquisition.FindAutofocus(SWAF_exp, timeoutSeconds=timeout)
152     except ApplicationException as e:
153         print('Application Exception : ', e.Message)
154     except TimeoutException as e:
155         print(e.Message)
156
157     print('Z-Position after initial SWAF : ', zSWAF)
158
159     return zSWAF

```

```

160
161 def checktableentry(datatable, entry2check='ImageSceneContainerName'):
162
163     num_col = datatable.ColumnCount
164     entry_exists = False
165     column = None
166
167     for c in range(0, num_col):
168         # get the current column name
169         colname = datatable.Columns[c].ColumnName
170         if colname == entry2check:
171             column = c
172             entry_exists = True
173             break
174
175     return entry_exists, column
176
177
178 def run_postprocessing(image, parameters={}, func='topography'):
179
180     if func == 'topography':
181
182         noise_low = parameters['noise_low']
183         noise_high = parameters['noise_high']
184         outputfolder = parameters['outputfolder']
185         ext = parameters['extension']
186
187         # converting to topo, with defined FZ noisecut
188         # in filter settings: 0-255 means no filter, 1-254 means cut one gray scale from top and from bottom
189         imgtop = Zen.Processing.Transformation.Topography.CreateTopography(image, noise_low, noise_high)
190         # saving file to the directory
191         topo_filepath = Path.Combine(outputfolder, Path.GetFileNameWithoutExtension(image.FileName) + ext)
192         Zen.Processing.Utilities.ExportHeightmapFromTopography(imgtop, topo_filepath)
193         print('Exported to : ', topo_filepath)
194         imgtop.Close()
195
196     return image
197
198 #####
199
200 # clear console output
201 Zen.Application.MacroEditor.ClearMessages()
202
203 # check the location of experiment setups and image analysis settings are stored
204 docfolder = Zen.Application.Environment.GetFolderPath(ZenSpecialFolder.UserDocuments)
205 imgfolder = Zen.Application.Environment.GetFolderPath(ZenSpecialFolder.ImageAutoSave)
206 imgfolder = Path.Combine(imgfolder, 'Guided_Acquisition')
207 # imgfolder = r'c:\Output\Guided_Acquisition'
208 format = '%Y-%m-%d_%H-%M-%S'
209
210 # get list with all existing experiments and image analysis setup and a short version of that list
211 expfiles = Directory.GetFiles(Path.Combine(docfolder, 'Experiment Setups'), '*.czexp')
212 ipfiles = Directory.GetFiles(Path.Combine(docfolder, 'Image Analysis Settings'), '*.czias')
213 expfiles_short = getshortfiles(expfiles)
214 ipfiles_short = getshortfiles(ipfiles)
215
216 # Initialize Dialog
217 GuidedAcqDialog = ZenWindow()
218 GuidedAcqDialog.Initialize('Guided Acquisition - Version : ' + str(version))
219 # add components to dialog
220 GuidedAcqDialog.AddLabel('1) Select Overview Experiment -----')
221 GuidedAcqDialog.AddDropDown('overview_exp', 'Overview Scan Experiment', expfiles_short, 0)
222 GuidedAcqDialog.AddCheckbox('fs_before_overview', 'OPTION - FindSurface (DF only) before Overview', False)
223 GuidedAcqDialog.AddCheckbox('SWAF_before_overview', 'OPTION - SWAF before Overview', False)
224 GuidedAcqDialog.AddIntegerRange('SWAF_ov_initial_range', 'Initial SWAF Range before Overview [micron]', 200, 50, 3000)
225 GuidedAcqDialog.AddLabel('2) Select Image Analysis to detect objects -----')
226 GuidedAcqDialog.AddDropDown('ip_pipe', 'Image Analysis Pipeline', ipfiles_short, 0)
227 GuidedAcqDialog.AddLabel('3) Select DetailScan Experiment -----')
228 GuidedAcqDialog.AddDropDown('detailed_exp', 'Detailed Scan Experiment', expfiles_short, 1)
229 GuidedAcqDialog.AddCheckbox('fs_before_detail', 'OPTION - FindSurface (DF only) before Detail', False)
230 GuidedAcqDialog.AddCheckbox('SWAF_before_detail', 'OPTION - SWAF before Detail', False)
231 GuidedAcqDialog.AddIntegerRange('SWAF_detail_initial_range', 'Initial SWAF Range before Detail [micron]', 100, 10, 1000)
232 GuidedAcqDialog.AddCheckbox('recallfocus_beforeDT', 'OPTION - Use RecallFocus (DF only) before Detail', False)
233 GuidedAcqDialog.AddLabel('4) Specify basefolder to save the images -----')
234 GuidedAcqDialog.AddFolderBrowser('outfolder', 'Basefolder for Images and Data Tables', imgfolder)
235
236 # show the window
237 result = GuidedAcqDialog.Show()
238 if result.HasCanceled:
239     message = 'Macro was canceled by user.'
240     print(message)
241     raise SystemExit
242
243 # get the values and store them
244 OverViewExpName = str(result.GetValue('overview_exp'))
245 ImageAS = str(result.GetValue('ip_pipe'))
246 DetailExpName = str(result.GetValue('detailed_exp'))
247 OutputFolder = str(result.GetValue('outfolder'))
248

```

```

249 fs_beforeOV = result.GetValue('fs_before_overview')
250 SWAF_beforeOV = result.GetValue('SWAF_before_overview')
251 SWAF_beforeOV_range = result.GetValue('SWAF_ov_initial_range')
252 fs_beforeDT = result.GetValue('fs_before_detail')
253 SWAF_beforeDT = result.GetValue('SWAF_before_detail')
254 SWAF_beforeDT_range = result.GetValue('SWAF_detail_initial_range')
255 RecallFocus = result.GetValue('recallfocus_beforeDT')
256
257
258 # print values
259 print('Overview Scan Experiment : ' + OverViewExpName)
260 print('Image Analysis Pipeline : ' + ImageAS)
261 print('Detailed Scan Experiment : ' + DetailExpName)
262 print('Output Folder for Data : ' + OutputFolder)
263 print('\n')
264
265 # check directory
266 OutputFolder = dircheck(OutputFolder)
267
268 # create a duplicate of the OVScan experiment to work with
269 OVScan_reloaded = cloneexp(OverViewExpName)
270
271 # active the temporary experiment to trigger its validation
272 OVScan_reloaded.SetActive()
273 time.sleep(hwdelay)
274 # check if the experiment contains tile regions
275 OVScanIsTileExp = OVScan_reloaded.IsTilesExperiment(blockindex)
276
277 ##### START OVERVIEW SCAN EXPERIMENT #####
278
279 if fs_beforeOV:
280     # initial focussing via FindSurface to assure a good starting position
281     Zen.Acquisition.FindSurface()
282     print('Z-Position after FindSurface: ', Zen.Devices.Focus.ActualPosition)
283
284 if SWAF_beforeOV:
285     zSWAF = runSWAF_special(OVScan_reloaded,
286                             delay=hwdelay,
287                             searchStrategy='Full',
288                             sampling=ZenSoftwareAutofocusSampling.Coarse,
289                             relativeRangeIsAutomatic=False,
290                             relativeRangeSize=SWAF_beforeOV_range,
291                             timeout=1)
292
293 # get the resulting z-position
294 znew = Zen.Devices.Focus.ActualPosition
295
296 # adapt the Overview Scan Tile Experiment with new Z-Position
297 if OVScanIsTileExp:
298     OVScan_reloaded.ModifyTileRegionsZ(blockindex, znew)
299     print('Adapted Z-Position of Tile OverView. New Z = ', '%.2f' % znew)
300
301 # execute the experiment
302 print('\nRunning OverviewScan Experiment.\n')
303 output_OVScan = Zen.Acquisition.Execute(OVScan_reloaded)
304 # For testing purposes - Load overview scan image automatically instead of executing the "real" experiment
305 # output_OVScan = Zen.Application.LoadImage(r'c:\Temp\input\OverViewScan_8Brains.czi', False)
306
307 # show the overview scan inside the document area
308 Zen.Application.Documents.Add(output_OVScan)
309 ovdock = Zen.Application.Documents.GetByName(output_OVScan.Name)
310
311 # save the overview scan image inside the select folder
312 output_OVScan.Save(Path.Combine(OutputFolder, ovsname))
313
314 ##### END OVERVIEW SCAN EXPERIMENT #####
315
316 # Load analysis setting created by the wizard or an separate macro
317 ias = ZenImageAnalysisSetting()
318
319 # for simulation use: 000 - RareEventExample.czis
320 ias.Load(ImageAS)
321
322 # Analyse the image
323 Zen.Analyzing.Analyze(output_OVScan, ias)
324
325 # Create Zen table with results for all detected objects (parent class)
326 AllObj = Zen.Analyzing.CreateRegionsTable(output_OVScan)
327
328 # Create Zen table with results for each single object
329 SingleObj = Zen.Analyzing.CreateRegionTable(output_OVScan)
330
331 # check for existence of required column names inside table
332 soi = SingleObj.GetBoundsColumnInfoFromImageAnalysis(True)
333
334 # 1st item is a bool indicating if all required columns could be found
335 columnsOK = soi.AreRequiredColumnsAvailable
336
337 if not columnsOK:

```

```

338     print('Execution stopped. Required Columns are missing.')
339     raise Exception('Execution stopped. Required Columns are missing.')
340
341 # show and save data tables to the specified folder
342 Zen.Application.Documents.Add(A11Obj)
343 Zen.Application.Documents.Add(SingleObj)
344 A11Obj.Save(Path.Combine(OutputFolder, 'OverviewTable.csv'))
345 SingleObj.Save(Path.Combine(OutputFolder, 'SingleObjectsTable.csv'))
346
347 # check the number of detected objects = rows inside image analysis table
348 num_POI = SingleObj.RowCount
349
350 ##### Prepare DetailScan #####
351
352 print('Starting DetailScan ...')
353
354 # create an duplicate of the DetailScan experiment to work with
355 DetailScan_reloaded = cloneexp(DetailExpName)
356
357 # active the temporary experiment to trigger its validation
358 DetailScan_reloaded.SetActive()
359 time.sleep(hwdelay)
360 # check if the experiment contains tile regions
361 DetailIsTileExp = DetailScan_reloaded.IsTilesExperiment(blockindex)
362
363 # test snap to change to the valid settings, e.g. the objective from the DetailScan
364 testsnap = Zen.Acquisition.AcquireImage(DetailScan_reloaded)
365 print('Acquire Test Snap using setting from DetailScan')
366 testsnap.Close()
367 # wait for moving hardware due to settings
368 time.sleep(hwdelay)
369
370 # move to 1st detected object
371 xpos_1st = SingleObj.GetValue(0, soi.CenterXColumnIndex)
372 ypos_1st = SingleObj.GetValue(0, soi.CenterYColumnIndex)
373 Zen.Devices.Stage.MoveTo(xpos_1st + dx_detector, ypos_1st + dy_detector)
374
375 if fs_beforeDT:
376     try:
377         # initial focussing via FindSurface to assure a good starting position
378         Zen.Acquisition.FindSurface()
379         print('Z-Position after FindSurface: ', Zen.Devices.Focus.ActualPosition)
380     except ApplicationException as e:
381         print('Application Exception : ', e.Message)
382         print('FindSurface (Definite Focus) failed.')
383
384 if SWAF_beforeDT:
385     zSWAF = runSWAF_special(DetailScan_reloaded,
386                             delay=hwdelay,
387                             searchStrategy='Full',
388                             sampling=ZenSoftwareAutofocusSampling.Coarse,
389                             relativeRangeIsAutomatic=False,
390                             relativeRangeSize=SWAF_beforeDT_range,
391                             timeout=0)
392
393 userecallfocus = False
394
395 if RecallFocus:
396     try:
397         # store current focus position inside DF to use it with RecallFocus
398         Zen.Acquisition.StoreFocus()
399         userecallfocus = True
400     except ApplicationException as e:
401         print('Application Exception : ', e.Message)
402         print('StoreFocus (Definite Focus) failed.')
403         userecallfocus = False
404
405
406 ##### START DETAILED SCAN EXPERIMENT #####
407
408 # get the actual Focus position
409 zpos = Zen.Devices.Focus.ActualPosition
410
411 # check for the column 'ID' which is required
412 ID_exists, column_ID = checktableentry(SingleObj,
413                                         entry2check='ID')
414
415
416 # execute detailed experiment at the position of every detected object
417 for i in range(0, num_POI, 1):
418
419     # get the object information from the position table
420     # get the ID of the object - IDs start with 2 !!!
421     #POI_ID = SingleObj.GetValue(i, 0)
422     POI_ID = SingleObj.GetValue(i, column_ID)
423
424     # get XY-stage position from table
425     xpos = SingleObj.GetValue(i, soi.CenterXColumnIndex)
426     ypos = SingleObj.GetValue(i, soi.CenterYColumnIndex)

```

```

427
428     # move to the current position
429     Zen.Devices.Stage.MoveTo(xpos + dx_detector, ypos + dy_detector)
430     print('Moving Stage to Object ID:', POI_ID, ' at :', '%.2f' % xpos, '%.2f' % ypos)
431
432     # try to apply RecallFocus (DF only) when this option is used
433     if userrecallfocus:
434         try:
435             # apply RecallFocus for the current position
436             Zen.Acquisition.RecallFocus()
437             zpos = Zen.Devices.Focus.ActualPosition
438             print('Recall Focus (Definite Focus) applied.')
439             print('Updated Z-Position: ', zpos)
440         except ApplicationException as e:
441             print('Application Exception : ', e.Message)
442             print('RecallFocus (Definite Focus) failed.')
443
444     print('New Z-Position before Detail Experiment will start:', zpos)
445
446     # if DetailScan is a Tile Experiment
447     if DetailIsTileExp:
448
449         print('Detailed Experiment contains TileRegions.')
450         # Modify tile center position - get bounding rectangle width and height in microns
451         bcwidth = SingleObj.GetValue(i, soi.WidthColumnIndex)
452         bcheight = SingleObj.GetValue(i, soi.HeightColumnIndex)
453         print('Width and Height : ', '%.2f' % bcwidth, '%.2f' % bcheight)
454         print('Modifying Tile Properties XYZ Position and width and height.')
455
456         # Modify the XYZ position and size of the TileRegion on-the-fly
457         print('Starting Z-Scan for current Object: ', '%.2f' % zpos)
458         print('New Tile Properties: ', '%.2f' % xpos, '%.2f' % ypos, '%.2f' % zpos, '%.2f' % bcwidth, '%.2f' % bcheight)
459         DetailScan_reloaded.ClearTileRegionsAndPositions(blockindex)
460
461         try:
462             DetailScan_reloaded.AddRectangleTileRegion(blockindex, xpos, ypos, bcwidth, bcheight, zpos)
463         except ApplicationException as e:
464             print('Application Exception : ', e.Message)
465
466     if not DetailIsTileExp:
467         print('Detailed Experiment does not contains TileRegions. Nothing to modify.')
468
469     # execute the DetailScan experiment
470     print('Running Detail Scan Experiment at new XYZ position.')
471     try:
472         output_detailscan = Zen.Acquisition.Execute(DetailScan_reloaded)
473     except ApplicationException as e:
474         print('Application Exception : ', e.Message)
475
476     # get the image data name
477     dtscan_name = output_detailscan.Name
478
479     """
480     Modification for multiscene images: container name needs to be part
481     of the filename of the detailed Scan.
482     First check if Column "Image Scene Container Name" is available
483     and then add Container Name to filename.
484     """
485
486     wellid_exist, column_wellid = checktableentry(SingleObj,
487                                                 entry2check='ImageSceneContainerName')
488
489     # in case Image Scene Container Name is not defined
490     if not wellid_exist:
491         message = 'Missing column ImageSceneContainerName in Image Analysis Results.\nPlease Select Features inside the Image
492         Analysis when needed.'
493         print(message)
494         container_name = 'empty'
495
496     # save the image data to the selected folder and close the image
497     output_detailscan.Save(Path.Combine(OutputFolder, output_detailscan.Name))
498     output_detailscan.Close()
499
500     # rename the CZI regarding to the object ID - Attention - IDs start with 2 !!!
501     if not wellid_exist:
502         # no wellID found inside table
503         newname_dtscan = 'DTScan_ID_' + str(POI_ID) + '.czi'
504     if wellid_exist:
505         # wellID was found inside table
506         well_id = SingleObj.GetValue(i, column_wellid)
507         newname_dtscan = 'DTScan_Well_' + str(well_id) + '_ID_' + str(POI_ID) + '.czi'
508
509     print('Renaming File: ' + dtscan_name + ' to: ' + newname_dtscan + '\n')
510     File.Move(Path.Combine(OutputFolder, dtscan_name), Path.Combine(OutputFolder, newname_dtscan))
511
512     ##### OPTIONAL POSTPROCESSING #####
513
514     if do_postprocess:
515         # do the postprocessing

```

```
515     image2process = Zen.Application.LoadImage(newname_dtscan, False)
516
517     # define the parameters for processing: Topography Export
518     parameters = {}
519     parameters['noise_low'] = 1
520     parameters['noise_high'] = 254
521     parameters['outputfolder'] = OutputFolder
522     parameters['extension'] = '.sur'
523
524     # run the processing and export and close the image
525     image2process = run_postprocessing(image2process, parameters=parameters, func='topography')
526     image2process.Close()
527
528 ##### END DETAILED SCAN EXPERIMENT #####
529
530     # restore the original OVScan experiment
531     OVScan_orig = Zen.Acquisition.Experiments.GetByName(OverViewExpName)
532     OVScan_orig.SetActive()
533
534     # delete the temporary experiments when all loops are finished
535     Zen.Acquisition.Experiments.Delete(OVScan_reloaded)
536     Zen.Acquisition.Experiments.Delete(DetailScan_reloaded)
537
538     # show the overview scan document again at the end
539     Zen.Application.Documents.ActiveDocument = ovdoc
540     print('All Positions done. Guided Acquisition Workflow finished.')
```

8 ToDo's and Limitations

The current version of this tool has still limitations and room for improvement.

- No yet built-in check to avoid double detailed scans for objects close to each other within one field of view.
- Currently only modification of rectangular tile regions are supported.

9 Disclaimer

This is an application note that is free to use for everybody. Use it on your own risk.

Especially be aware of the fact that automated stage movements might damage hardware if the system is not setup properly and not correctly calibrated.

Please check everything in simulation mode first!

Carl Zeiss Microscopy GmbH's ZEN software allows connection to the third party software, Python. Therefore Carl Zeiss Microscopy GmbH undertakes no warranty concerning Python, makes no representation that Python will work on your hardware, and will not be liable for any damages caused by the use of this extension.

By running and using this script confirm that you understood the involved risks you agree to this disclaimer.