


# ZEN and APEER - Open Ecosystem for integrated Machine-Learning Workflows




**Sebastian Rhode**

Product Owner Machine Learning, Technology Center Software, Munich




- **Deep Learning** is a special subset of general **Machine Learning**, which is again a subset of **Artificial Intelligence (AI)**.
- Machine learning gives computers the ability to **learn without being programmed** explicitly for a specific task
- Two main categories exist in Machine Learning: Supervised and Unsupervised methods

# Deep Learning – What are Artificial Neural Networks?



- Deep Learning itself consists of methods to analyze data using an approach that mimics the way a person would try to make sense of it
- Deep Learning methods use a layered structure of algorithms called artificial neural networks (ANN)
- ANNs have an input and output layer
- Layers in between are hidden layers → their values cannot be observed (easily) during the training
- ANN gets "deeper" if the number of hidden layers increases

# Machine Learning – Deep Learning versus Machine Learning (simplified)




- Machine Learning algorithms are rather "simple" compared to ANNs
- classical Machine Learning **methods are using features typically "engineered" by human intervention**
- Deep Learning **methods extract features automatically** and learn which features work best
- Deep Learning typically requires a lot more data due to the complex algorithms used

# Application of Machine Learning methods in Microscopy


## Image Analysis and Processing

### Classification


Classifying an entire image or individual objects



Slide vs 96 WellPlate (CD7 or SampleFinder)




Classify Objects (ZEN Object Classification)




Good part vs anomaly

### Segmentation

Refers to classifying at a pixel level



Semantic Segmentation



Instance Segmentation

### Processing

Image corrections and enhancements




Denoising




Super-resolution


# Machine-Learning - The complexity of the task defines the approach



Complexity of  
the application



Binary segmentation




Nuclear segmentation

# amount of  
data required

single image

several images  
( $\geq 2$  images)



Segmenting Membranes

# parameters  
for training

1

10

100

$1e+3$

$1e+4$

$1e+5$

$1e+6$


$1e+7$

**Thresholding**

**Pixel Classification**

**Deep Neural Networks**

# Machine Learning – Train a Model and use a Model



*“I can get decent images from anywhere but I am looking for information”*

- Potential customer during her visit to Pleasanton

## Segmentation and Analysis are today's biggest problems.

A powerful and reproducible segmentation is the mandatory precursor for any downstream image analysis and measurements.

The actual purpose of an image is “to be analyzed”.

- Images can be “hard” or challenging to segment because of their nature.
- Further sorting and classifying segmented objects can be difficult.
- Users are often no experts and need a simple tool.

# ZEN and APEER – Open Ecosystem for integrated Machine-Learning Workflows

Trained model from Segmentation from “anywhere”



apeer




Train a DNN model

Train a DNN model


Train a  
PixelClassifier  
in Intellesis

Trained Model

Get a trained  
DNN model as  
a service



# ZEISS – Open Ecosystem for integrated Machine-Learning workflows



- APEER and ZEN are part of the same ecosystem
  - Benefit from the Annotation and Training tools
  - Use trained neural networks inside APEER workflows
- Creating a model on can be done **by users** or as a **paid service** (SW-SCS)
- **The “exchange” currency inside this ecosystem is the trained model**
- ML models developed on APEER (or anywhere) can be used in integrated workflows
  - Use inside ZEN workflows like Image Analysis, Guided Acquisition, Material Module etc.
  - Benefit for support for large data sets, Tiling Capabilities and Integration

# ZEN Intellesis – Open Ecosystem for integrated Machine-Learning Workflows

Use your own TF2 models in ZEN – Open Source PyPi package CZMODEL



**czmodel 1.1.0**

`pip install czmodel`

Released: Mar 8, 2021

A conversion tool for TensorFlow ANNs to CZModel

**Navigation**

- Project description** (selected)
- Release history
- Download files

**Project description**

This project provides simple-to-use conversion tools to generate a CZModel file from a [TensorFlow](#) model that resides in memory or on disk to be usable in the [ZEN Intellesis](#) module starting with ZEN blue >=3.2 and ZEN Core >3.0.

This version of czmodel produces the following model version: 3.1.0

Please check the following compatibility matrix for ZEN Blue/Core and the respective version of a czmodel file. Version compatibility is defined via the [Semantic Versioning Specification \(SemVer\)](#).

Model	ZEN Blue	ZEN Core
3.1.0	> 3.3	> 3.2
3.0.0	> 3.1	> 3.0

If you encounter a version mismatch when importing a model into ZEN, please check for the correct version of this package.

**System setup**

The current version of this toolbox only requires a fresh Python 3.x installation. It was tested with Python 3.7 on Windows.


**Meta**

License: Apache Software License


- Open and standardized “container” to store ML models and PyPi
- The CZMODEL is the “glue” or “exchange format” between SW tools for Machine-Learning at ZEISS
- No new model format!
- Support for TF2.SavedModels and ONNX models
- Internally used by ZEN blue, ZEN core, APEER-ML and other parts of ZEISS (even beyond Microscopy)
- Free to use for everybody
- Allows external data scientists to integrate their own models into ZEN workflows

# ZEN Intellesis – Open Ecosystem for integrated Machine-Learning Workflows

Use your own TF2 models in ZEN – Open Source PyPi package CZMODEL



# ZEN Intellesis – Open Ecosystem for integrated Machine-Learning Workflows



- **Simple User Interface for Labelling and Training**

- Label your datasets using a clean and simple UI, train a Pixel Classifier (see: [Feature Extractors](#))

- **Integration into ZEN Measurement and Processing Framework**

- Use trained models easily inside complete Image Analysis pipelines, Material Modules or inside scripted workflows

- **Open Platform – Import your own trained networks**

- Import your own model and use it seamlessly integrated in ZEN workflows and benefit from Image Tiling & Fusion strategy (see: [Importing external networks into ZEN](#))

- **Support for Multi-dimensional Datasets**

- Import any multi-dimensional dataset incl. 3<sup>rd</sup> party file formats from other vendors



apeer




ONNX

K Keras



scikit-learn

# Simple User Interface for Labeling and Training



- Simple and intuitive User Interface to train
  - Segmentation models (APEER and ZEN)
  - Object Classifier (ZEN)
- ZEN Software with Image Analysis Tool, Scripting Interface and pre-defined workflows
- AppManager Infrastructure and ZeissPy distribution integrated in installer
- Tiling Strategy and Image Fusion for nD datasets incl. support for 3rd party data formats
- Open-Source Package `czmodel`: <https://pypi.org/project/czmodel/>
- Import of externally trained deep-learning models into ZEN platform
- **Support for TensorFlow2 and ONNX models**
- Cloud-based Annotation and Training Platform for neural networks: [www.apeer.com](http://www.apeer.com)
- **Data Annotation and Model Training as a Service**



- **Multispectral Feature Extraction** – all channels will be used to segment a pixel
- **Class Segmentation** – hierarchical structures with independent segmentation per object class
- **Engineered Feature Sets and Deep Feature Extraction and pre-trained networks**
- **Random Forrest Classifier + DNN via TF.SaveModels**
- **Public [DNN Model Specification](#) on [ZEISS Microscopy GitHub](#) page**
- Post processing: **Conditional Random Fields and confidence thresholds**
- IP-Functions for **Creating Masks** and **Scripting Integration**
- Client-Server Architecture (REST-API) with **Client-side Tiling Strategy**

```
def runmodel(image, model,
            use_confidence=False,
            confidence_threshold=0,
            format='MultiChannel',
            extractclass=False,
            class2extract_id=0,
            addseg=False,
            adapt_pixeltypes=True):

    # define the desired output format
    if format == 'MultiChannel':
        # output image will have one channel per model class
        segf = ZenSegmentationFormat.MultiChannel
    if format == 'Labels':
        # output image will have one channel with distinct labels per model class
        segf = ZenSegmentationFormat.Labels


    # classify pixels using a trained model incl import deep-learning models
    if use_confidence:
        try:
            # run the segmentation and apply confidence threshold to segmented image
            outputs = Zen.Processing.Segmentation.TrainableSegmentationWithProbabilityMap(image, model, segf)
            seg_image = outputs[0]
            conf_map = outputs[1]
            print('Apply Confidence Threshold to segmented image.')
            seg_image = Zen.Processing.Segmentation.MinimumConfidence(seg_image, conf_map, confidence_threshold)
            conf_map.close()
            del outputs
        except ApplicationException as e:
            seg_image = None
            print('Application Exception : ', e.Message)

    if not use_confidence:
        try:
            # run just the segmentation
            seg_image = Zen.Processing.Segmentation.TrainableSegmentation(image, model, segf)
        except ApplicationException as e:
            seg_image = None
            print('Application Exception : ', e.Message)


    return seg_image
```

- **Labels is valid for all channels of an image** and the information from all channels will be used to classify that pixel (X,Y)
- **Label as few pixels as possible** and rather do more iterations of training.
- Sometimes it is helpful to **use more classes in order to get better results for each class**.
- **Labeling large homogenous areas with many pixels will slow down the training** but give not necessarily give better results.
- It is **good practice to label “roughly” the same number of pixels for every class**, e.g., try to avoid labeling much, much more pixels for one class compared to another class
- The **default feature extractors are using the fewest number of parameters** (start here first) while the Deep Feature extractors use the largest number. Use those only when the Default one do not work well.
- If those methods do not work, **consider training a Deep Neural Network on APEER or consider getting a model as a service**.



# Intellesis Segmentation - Workflows and Measurement Integration




# Actionable Information – Where can trained models be used




# Training UI in ZEN blue




# Integration into ZEN Image Analysis and Software Modules



# Integration into Image Analysis



# Integration into Image Analysis



- Intellesis Segmentation allows to use 7 different feature extractors to create the feature vector that will be used for the subsequent Pixel Classification using Random Forrest.
- For calculating the features of one specific pixel, various filters with various filter sizes and parameters are applied to the region around this pixel.
- All the filter results are concatenated yielding the final feature vector describing the pixel. Considering the filter sizes of all filters we obtain a dimensionality of 25/33 for the Basic Features 25/33
- For Deep Feature Extraction ZEN Intellesis uses the results of an intermediate layer of the VGG19 network as input features for the Pixel Classification using a Random Forest Classifier.

### For more detailed Information:

#### [Intellesis Segmentation - Feature Extraction for Pixel Classification](#)

Remark: Do not confuse Intellesis Segmentation via PixelClassification with using pre-trained deep neural networks, which is also possible with Intellesis Segmentation.


- Each decision tree in the forest considers a subset of features when forming “questions” and only has access to a subset set of the training data
- This increases diversity in the forest leading to more robust overall predictions and the name ‘**Random Forest.**’
- Prediction = average of all the individual decision tree estimates
- Regression task → predicting a continuous value of ...
- **Classification task → discrete class labels such as A or B**
  - Random Forest will take a majority vote for the predicted class

<https://medium.com/@williamkoehrsen/random-forest-simple-explanation-377895a60d2d>


[scikit-learn - Random Forest](#)

[Wikipedia - Random Forrest](#)

# Random Forest Classifier – Decision Trees




## Postprocessing – Conditional Random Fields (CRF)




- CRF tries to create smooth and sharper edges by re-classifying pixels based on confidence levels in their neighborhood
- after CRF there is no real probability map from the Random Forrest anymore, all confidence values will be set to 100%
- For more details: [PostProcessing - Conditional Random Fields](#)

# ZEN Intellesis Object Classification – Training UI



# ZEN Intellesis Object Classification – Result Visualization



- **Simple User Interface for Labelling Objects and Train a Classifier**

- Label the segmented objects using a clean and simple UI and train an Object Classifier using already analyzed images

- **Integration into ZEN Measurement and Processing Framework**

- Use trained Object Classifier inside complete Image Analysis pipelines, Material Modules or inside scripted workflows

- **Support for Multi-dimensional Datasets**

- Import any multi-dimensional dataset incl. 3<sup>rd</sup> party file formats from other vendors



apeer

ANACONDA

python™

TensorFlow

ONNX

K Keras

DASK

scikit-learn

- **Requires an analyzed image** with already segmented objects
  - segmentation method can be "anything"
  - Thresholds, PixelClassifier, DeepNeuralNetwork
- **Multispectral Feature Extraction** – all channels can be used to extract intensity-based features
- **Geometrical and Shape Features** – all available features from the Image Analysis will be used automatically
- **Random Forrest Classifier** to process the feature table
- IP-Functions for Classification and **Scripting Integration for Automation**

```
# get all available object classification models
ocmodels = ZenIntellesis.ObjectClassification.ListAvailableModels()

for obcmode in ocmodels:
    print('TOC Model Name:', obcmode.Name)
    print('Description   :', obcmode.Description)
    print('Status        :', obcmode.Status)

    # get training images
    training_images = obcmode.TrainingImages
    for ti in training_images:
        print(ti)

# select a specific model
myobjmodel = find_objclassmodel('my_toc_model')
print('Use model:', myobjmodel.Name)

# import a model
model2import = r"C:\Object_Classification_Models\model1.cztoc"

# classify an analyzed image (inplace)
myobjmodel.Classify(image, appendFeatures=False) # Where is append feature option

# define IA class of interest
IAclass_all = 'Statistics objectA'
IAclass_single = 'objectA'

# define the object class of interest
objclass = 'alive'


# export table for ALL objects (statistics)
table_all = Zen.Analyzing.CreateRegionsTable(image, IAclass_all)
Zen.Application.Documents.Add(table_all)

# export table for SINGLE objects (the ones classified by the model)
table_single = Zen.Analyzing.CreateRegionTable(image, IAclass_single)


# extract a specific object class for from table by filtering (in place)
table_single.RemoveRows(lambda row: row[myobjmodel.Name] != objclass)
```

- ZEN Intellesis Object Classification is a **dedicated SW module** for object classification, and it is independent from ZEN Intellesis Segmentation.
- Its task is to **classify (already segmented) objects based on measured features**, not to measure parameters.
- Object Classification uses **intensity-based features** and all available **geometrical and shape-based features**.
- In other words, **object classification “sorts” objects** of one specific region class from image analysis (cells) into different type of objects (cell type A and type B).
- It **requires analyzed images (!) for training and for classification** and can be used after ZEN Image Analysis (not within) – it a subsequent workflow step.
- The **type of segmentation does not matter** – any methods (classical or machine-learning) can be used to create the segmented objects.
- The **input** for the classifier is the **table with measured features**, not the actual image itself.

# Intellesis Object Classification - Workflows and Measurement Integration



# ZEN Intellesis Segmentation and Classification Workflow



# ZEN Machine-Learning – Full integration into Python Scripting Interface



```
def runmodel(image, model,
             use_confidence=False,
             confidence_threshold=0,
             format='MultiChannel',
             extractclass=False,
             class2extract_id=0,
             addseg=False,
             adapt_pixeltype=True):

    # define the desired output format
    if format == 'MultiChannel':
        # output image will have one channel per model class
        segf = ZenSegmentationFormat.MultiChannel
    if format == 'Labels':
        # output image will have one channel with distinct labels per model class
        segf = ZenSegmentationFormat.Labels

    # classify pixels using a trained model incl import deep-learning models
    if use_confidence:
        try:
            # run the segmentation and apply confidence threshold to segmented image
            outputs = Zen.Processing.Segmentation.TrainableSegmentationWithProbabilityMap(image, model, segf)
            seg_image = outputs[0]
            conf_map = outputs[1]
            print('Apply Confidence Threshold to segmented image.')
            seg_image = Zen.Processing.Segmentation.MinimumConfidence(seg_image, conf_map, confidence_threshold)
            conf_map.Close()
            del outputs
        except ApplicationException as e:
            seg_image = None
            print('Application Exception : ', e.Message)

    if not use_confidence:
        try:
            # run just the segmentation
            seg_image = Zen.Processing.Segmentation.TrainableSegmentation(image, model, segf)
        except ApplicationException as e:
            seg_image = None
            print('Application Exception : ', e.Message)

    return seg_image
```

```
# get all available object classification models
ocmodels = ZenIntellessis.ObjectClassification.ListAvailableModels()

for obcmode in ocmodels:
    print('TOC Model Name:', obcmode.Name)
    print('Description   :', obcmode.Description)
    print('Status       :', obcmode.Status)

    # get training images
    training_images = obcmode.TrainingImages
    for ti in training_images:
        print(ti)

# select a specific model
myobjmodel = find_objclassmodel('my_toc_model')
print('Use model:', myobjmodel.Name)

# import a model
model2import = r"C:\Object_Classification_Models\model1.czto"

# classify an analyzed image (inplace)
myobjmodel.Classify(image, appendFeatures=False) # Where is append feature option

# define IA class of interest
IAclass_all = 'Statistics objectA'
IAclass_single = 'objectA'

# define the object class of interest
objclass = 'alive'

# export table for ALL objects (statistics)
table_all = Zen.Analyzing.CreateRegionsTable(image, IAclass_all)
Zen.Application.Documents.Add(table_all)

# export table for SINGLE objects (the ones classified by the model)
table_single = Zen.Analyzing.CreateRegionTable(image, IAclass_single)




# extract a specific object class for from table by filtering (in place)
table_single.RemoveRows(lambda row: row[myobjmodel.Name] != objclass)
```

# Clearly defined steps to guide the user from annotation to training to segmentation

Annotate on APEER

Train your specific U-net

Segment other images –  
on APEER or in Intellesis



# Clearly defined steps to guide the user from annotation to training to segmentation




Overview of APEER ML

## Overview of APEER ML

Takes you from image annotation to training and segmentation using machine learning tools.

### 3 steps to image segmentation

The journey from images to segmentation goes through 3 primary steps, annotate, train and segment. Follow the full process or jump to a specific step.



An image showing a microscopy slide with various cellular structures. A small circular play button is overlaid on the image.

**Annotate**  
Upload your images, annotate and save labels to your APEER storage. Remember that you have 100GB free storage on APEER.

[Get started](#)



An image showing a computer monitor displaying a software interface with multiple panels and data visualizations. A small circular play button is overlaid on the image.

**Train**  
Use your images and saved labels to train a model for semantic (pixel) or instance (object) segmentation.

[Get started](#)



An image showing a microscopy slide with segmented regions highlighted in different colors (blue, green, red). A small circular play button is overlaid on the image.

**Segment**  
Use your trained model to automate the image segmentation process. APEER takes care of image handling and hyperparameter tuning so you can focus on the results.



[Get started](#)

**Gain insights, do not stop with segmentation!**

The 'Segment' step includes semantic/instance segmentation followed by object measurements. The result contains a comprehensive collection of parameters, including centroid, diameter, and moments. For additional tools, please browse through our [public modules](#).






# Clearly defined steps to guide the user from annotation to training to segmentation



# Clearly defined steps to guide the user from annotation to training to segmentation

User will never see nor must interact with these complex background tasks.







## Flexible & Hierarchical Segmentation Scheme (Class Segmentation)


- Cell - classical threshold
  - Nucleus – DNN Nucleus Detector
    - Object A - trained model A
    - Object B - trained model B
  - Object C - classical threshold

# Class Segmentation – Hierarchical structures with independent segmenter per class


- Import dialog for channel selection incase of multichannel image
- Single Channel is the default
- Next import allows only selecting the channel



# Class Segmentation – Hierarchical structures with independent segmenter per class




# Class Segmentation – Hierarchical structures with independent segmenter per class



# Nucleus Detector – UNet architecture (vgg16 encoder)

Try out pretrained models from ZEISS



BGR to RGB



Per-image-standardization



Input



Conv 2D



Max-Pooling 2 x 2



Upsampling



Concatenate



Batch Normalization





Output (Conv2D + Softmax)

- Number above represent the array sizes
- numbers below represent the number of feature maps (per layer)


# Use Deep Neural Networks for a specific application



- Import pretrained networks provided ZEISS, APEER, SCS-Team or Researchers
- ZEISS example application: Pre-trained Nucleus Detection
- Use directly for Segmentation, Image Analysis incl. Class Segmentation



# Pretrained Network – Integration into ZEN Image Analysis




- Network can be used to create a fully functional Image Analysis Pipeline by one single click
- Such IA settings will now also have the advanced binary postprocessing options
- Confidence Threshold is available

- Starting with ZEN Blue 3.1 it is possible to download / import pretrained networks for Segmentation
- Check it out: [Model Download](#)
- While Thresholding combined with classical PreProcessing function is really fast, it can be challenging (or impossible), especially for the non-expert user, to get the same results
- To illustrate this the following slides shows the same image analyzed by an Image Analysis Setting using
  - preProcessing + Thresholding
  - DNN-based Segmentation


## Side Note – Deep Neural Network versus Threshold

Same image data – different display curve




# Side Note – Deep Neural Network versus Threshold

Same image data – different results




# Advanced Class Segmentation incl. trained models

Plugin your model into any image analysis class





# Advanced Class Segmentation incl. trained models

## Easily visualize the result of your machine-learning based segmentation



# Nucleus Detector – Use inside Class Segmentation








Sample courtesy of Dr. Michael Sodek, CRUK Institute, London

- XRM scan: Drosophila
- Elyra SIM: Neuron – spines & dendrite
- Celldiscoverer 7 (BF): Organoid and migrated cells
- EM: classification of cell compartments
- Celldiscoverer 7 (BF + FL): Scratch assay – area and mitotic cells
- Airyscan Fast: Muscle section – different myosin types


## Application Examples – Life Science





XRM scan: Drosophila




Neuron – spines & dendrite





Organoid and migrated cells




Classification of cell compartments






Scratch assay – area and mitotic cells




Muscle section – different myosin types



## Application Examples – Life Science



- Cells image using Phase-Gradient Contrast on a CD7
- Labeled with 2 classes inside Intellesis Training UI
- Feature Extractor:  
DeepFeatures256 + CRF Postprocessing



- Cells image using Phase-Gradient Contrast on a CD7
- Nucleus Segmentation using a Deep Neural networks trained on the APEER platform
- Segmented nuclei shown inside ZEN blue Image Analysis

[Intellessis Landing Page](#)

Get the trial license and try it out!

---

[APEER - Machine Learning Overview](#)

---

Additional Content on GitHub:

[Machine Learning](#)

[Machine Learning - Feature Extractors](#)

[Post Processing](#)



We make it visible.